

CELOXTM **RAPID GAUZE**

Stops Bleeding. Saves Lives.

CeloxTM Rapid The fastest acting battlefield haemostat

- Works with just 60 seconds compression¹
- Stops life-threatening bleeding - fast
- Rapid action reduces blood loss²



Contact us at www.celoxmedical.com




Celox™ Rapid - The new standard for haemorrhage control

Time is critical when treating severe injuries in the field and saving time increases the chances of survival for the casualty in an emergency or hostile situation. To be effective, haemostats need to work fast.

The current generation of haemostatic agent require at least 3 minutes of compression to work. Celox™ Rapid Gauze removes the delay by speeding up packing time and reducing compression time to just **60 seconds**, while using the proven haemostat technology of Celox™.

How Celox™ Rapid works


Celox™ Rapid contains Chito-R™ activated chitosan granules that stick to wet tissue in the wound³. Chito-R™ is bonded to a high density gauze that packs faster than traditional gauze.



- Chito-R forms an adhesive gel that seals the wound to stop the flow of blood, speeding up compression time;
- Celox™ Rapid absorbs blood fluid, swells and forms a gel;
- Works on hypothermic blood;
- No minerals, chitosan is made up of glucosamine and N-acetyl glucosamine.

Cut treatment time in half

Treatment time (minutes)



Because Celox™ Rapid sticks to wet tissue, it only needs 60 seconds of compression with no loss of performance; as opposed to three minutes for other haemostats.

What is more, Celox™ Rapid also saves time packing the wound⁴, so the real saving is more than two minutes.

What this means

- Celox™ Rapid has been independently shown to reduce blood loss by 60% compared to Quikclot Gauze*²
- Celox™ Rapid stays in place during evacuation⁵
- In an emergency, speed means safety for the whole team
- Faster treatment improves casualty survival and helps treatment of multiple casualties

Ordering information

Z-Fold Gauze: FG08839021
NSN: 6510 99 184 6453

Medtrade Products Ltd.

Electra House, Crewe Business Park, Crewe CW1 6GL UK
T: +44 (0)1270 500019 F: +44 (0)1270 500045

MT/13/003

www.celoxmedical.com

