

# EMERGENCY'S™

## Blue light bag EM13.014

Organized and handy transport for emergencies



**Practical and light emergency's bag with two adjustable inner separators and a rack pocket in the flap to ease the organization of the sanitary material. It also have three external pockets, two of them are diaphanous and one is in the front and is internally subdivided and with elastic rubber-bands.**

### Features

- High resistance zips
- Night & day reflecting strips
- Washable material
- Reinforced handles and shoulder straps

### Technical data

- Measure Main Body: 44 x 25 x 27 cm.
- Total Volume: 21,60 L
- Material: 600D Polyester with PVC .
- Weight: 1,29 Kg

# EMERGENCY'S™

websites

 ELITE BAGS

 MILITARY

 EMERGENCY'S™

tee-uu

social media

