# EMERGENCY'S"

## Light bag Tarpaulin

#### EM13.021

Light and compact emergency bag developed in Tarpaulin, designed to facilitate cleaning and ensure water, body fluids and fluid repellency in general. It has two adjustable internal dividers and a mesh pocket on the flap to facilitate the organization of the sanitary material, as well as three external pockets, two of them diaphanous and one frontal sub divided in its interior and with elastic bands.



#### Features

- High resistance zips.
- Night & day reflecting strips.
- Washable material.
- Reinforced handles and shoulder straps

### **Technical data**

Weight	1.23 Kg
Dimensions	44 x 25 x 27 cm
Material	TARPAULIN
Color	RED
Capacity (Liters)	21.60

Vial de los Cipreses, 9, 03006 Alicante, Spain. Telf: (+34)902 212 210 · info@elitebags.es · www.elitebags.es



