

EMERGENCY'S™

Great Capacity Bag

EM13.012

Basic Life Support padded and large capacity Bag, designed to house diverse sanitary material. Subdivided in its interior with 4 medium compartments and one large capacity elongated compartment. Transparent pocket in the inner flap, designed in order to easily see and organise different material, thanks to the elastic rubber-bands which guarantee a proper fastening. It also have five external pockets, two side pockets, one on the back and one on the top which are diaphanous, and finally a frontal pocket which can be internally divided thanks to pockets and rubber-bands.



Features

- Tarpaulin made base (waterproof material)
- High resistance zips
- Night & day reflecting strips
- Washable material
- Reinforced handles and shoulder straps
- External card holder
- External top transparent pocket to keep list of sanitary material.

Technical data

Weight	1.76 Kg
Dimensions	55 x 31.5 x 33 cm
Material	600 D POLIESTER WITH PVC
Color	BLUE
Capacity (Liters)	37.83