EMERGENCY'S™

Rescue Backpack

EM13.027

Medium capacity Basic Life Support (BLS) Backpack. Its opening with folding flap allows access to the tools necessary for the performance in different emergency situations, placing them into a single workspace.





Features

- Ergonomic main handle with velcro closing.
- Large pockets in both sides.
- · Adjustable padded shoulders and backpack.
- Padded back pannel.
- Large front pocket ready to accommodate the defibrillator.
- Interior drop for convenient and easy access inside instrumental.
- Internal compartments with mesh and transparent window.

Technical data

Weight 1.22 Kg

Dimensions 21 x 37 x 45 cm

Capacity (Liters) 29.75

Color ORANGE

Material 600 D POLIESTER WITH PVC







