



## MIL-TEC Υδροδοχείο Πλάτης με Ρυθμιζόμενους Ιμάντες Ώμου και Ζώνη Στήθους 3lt με D-rings - Άμμου

€32.18

[Αγορά](#)

### Περιγραφή

## MIL-TEC Υδροδοχείο Πλάτης με Ρυθμιζόμενους Ιμάντες Ώμου και Ζώνη Στήθους 3lt με D-rings - Άμμου

1000D Cordura Nylon (ΑΡΙΣΤΗ ΠΟΙΟΤΗΤΑ - ΜΕΓΙΣΤΗ ΑΝΤΟΧΗ)  
Περιλαμβάνεται δοχείο 3lt.

Το σύστημα ενυδάτωσης φέρει 6 D-rings π.χ. για να συνδέεται με σακίδιο.

- Σακίδιο ενυδάτωσης με ρυθμιζόμενους ιμάντες ώμου και ζώνη στήθους.
- Πόσιμο σωλήνα με βαλβίδα.
- Ιμάντας μεταφοράς.
- Μεγάλο άνοιγμα για εύκολο γέμισμα και άδειασμα.
- Ο σωλήνας πόσης του δοχείου είναι επικαλυμμένος με neopren

## Mil-Tec Water Pack Mil Spec Coyote

Quality hydration pack with water bladder and adjustable shoulder straps, equipped with 6 D-rings that allows connection with other rucksacks or attaching extra equipment if needed. It goes with leak proof and taste free water carrier with large opening for easy filling and emptying, a tube insulated with neoprene cover, convenient On/Off valve, hygienic interchangeable mouthpiece and a cap.

- Capacity: 3.0 litres

- Mesh pocket on the back
- Adjustable shoulder straps
- Adjustable waist strap
- Adjustable chest strap
- Grab handles on the top and front
- Drainage hole at the bottom
- D-rings
- Material: 1000D Cordura nylon
- Dimensions: approx. 23x46cm (9"x18")
- Weight: 605g
- Hydration bladder included (Microban)
- Manufacturer: Mil-Tec

---

The document was created on: Ιανουάριος 10, 2025