



## MIL-TEC Υπνόσακος "PILOT"

**€24.11**

[Αγορά](#)

## Περιγραφή

Upper: 100% Polyester

Lining: 100% Cotton

Filling: 100% Polyester

Transportation bag: 100% Polyester

Dimensions: 185 x 75 cm

Pack size: 40 x 24 x 29 cm

Weight: 1700 g

Extreme range: until ca. 0 degrees

- swen pack bag (also usable as head pad)
- sleeping bag also usable as a blanket
- 2-way-zipper
- zipper-stick-protection (hook and loop closure)
- warmth collar along the zipper