



## MIL-TEC Θερμός Κενού Αέρος - 0,5 lt & 1 lt

€7.66

[Αγορά](#)



## Περιγραφή

Ένα θερμό κενού αέρος από την εταιρία MIL-TEC. Κατασκευασμένο από υψηλής ποιότητας μέταλλο που το καθιστά ανθεκτικό σε χτυπήματα και αντίξοες συνθήκες.

Χωρητικότητα	Βάρος	Ύψος	Διάμετρος
0,5 λίτρα	345 gr	24,5 cm	7 cm
1 λίτρο	600 gr	32 cm	8cm

- Άθραυστο μπουκάλι
- Το καπάκι μπορεί να χρησιμοποιηθεί και σαν ποτηράκι
- Διαθέτει μια πρακτική βαλβίδα η οποία βοηθά στο να χύνεται ορισμένη ποσότητα από το μπουκάλι και έτσι να μπορεί να χρησιμοποιηθεί με το ένα χέρι.
- Μεγάλο άνοιγμα για εύκολο καθαρισμό
- Μπορεί να χρησιμοποιηθεί για κρύα αλλά και ζεστά ροφήματα.

## OD STAINLESS STEEL THERMO BOTTLE

High-grade steel vakuuum thermo bottle.

Total volume Weight Height Diameter

0,5 liters	345 gr	24,5 cm	7 cm
1 liters	500 gr	32 cm	8 cm

- unbreakable thermo bottle with flat outside
- drinking cup as sealable cap
- practical quick stop pouring out system for comfortable single-hand control
- wide opening for cleaning
- for hot and cold drinks

---

The document was created on: Ιανουάριος 9, 2025