



ΜΑΧΑΙΡΙ ΑΛΒΑΙΝΟΧ WOODEN, STAMINA, 32273

€39.11

[Αγορά](#)

Περιγραφή

ΜΑΧΑΙΡΙ ΑΛΒΑΙΝΟΧ WOODEN, STAMINA, 32273

Registered model
Designed in Spain
Handle: Decorated CNC stamina
Blade: Stainless steel
Blade size: 10.5 cm
Total size: 22.5 cm
Blade thickness: 4.8 mm
Bolster: Stainless steel
Leather sheath

Χαρακτηριστικά

The document was created on: Ιανουάριος 9, 2025