



ΜΑΧΑΙΡΙ Κ25, TACTICO, KRAKEN, Titanium Coated, 32373

€19.76

[Αγορά](#)

Περιγραφή

ΜΑΧΑΙΡΙ Κ25, TACTICO, KRAKEN, Titanium Coated, 32373

Mod. Regist.©
Designed in Spain
Brand: K25
Model: KRAKEN
Tactical knife
Handle: G10
Blade: Stainless steel (Titanium coated)
Blade thickness: 3.3 mm
Blade size: 7.2 cm
Total size: 14 cm
Kydex sheath

Χαρακτηριστικά

The document was created on: Ιανουάριος 9, 2025