



ΜΑΧΑΙΡΙ Κ25, TACTICO, G10 Handle, 32371

€19.76

[Αγορά](#)

Περιγραφή

ΜΑΧΑΙΡΙ Κ25, TACTICO, G10 Handle, 32371

Mod. Regist.©

Designed in Spain

Brand: K25

Tactical knife

Handle: G10

Blade: Stainless steel (Titanium coated)

Blade thickness: 3 mm

Blade size: 7.5 cm

Total size: 14 cm

Kydex sheath

Χαρακτηριστικά

The document was created on: Ιανουάριος 9, 2025