



ΜΑΧΑΙΡΙ Κ25, Tactical Knife, Titanium Coated ,MI-28 23.5cm, 32391

€34.68

[Αγορά](#)

Περιγραφή

ΜΑΧΑΙΡΙ Κ25, Tactical Knife, Titanium Coated ,MI-28 23.5cm, 32391

Designed in Spain

Brand: K25

Model: MI-28

Tactical knife

Handle: Rubber

Blade: Stainless steel (Titanium coated)

Blade thickness: 3.8 mm

Blade size: 11.5 cm

Total size: 23.5 cm

Rigid pouch

Χαρακτηριστικά

The document was created on: Ιανουάριος 9, 2025