



MAXAIPI K25, Tactical Knife, Titanium Coated ,Coyote Cord Wrapped, 32378

€31.05

[Αγορά](#)

Περιγραφή

MAXAIPI K25, Tactical Knife, Titanium Coated ,Coyote Cord Wrapped, 32378

Designed in Spain

Brand: K25

Handle: Cord wrapped + G10

Blade: Stainless steel (Titanium coated)

Blade thickness: 5.8 mm

Total size: 25 cm

Nylon sheath

Χαρακτηριστικά

The document was created on: Ιανουάριος 9, 2025