



ΜΑΧΑΙΡΙ Κ25, PYTHON, 32338

€19.27

[Αγορά](#)

Περιγραφή

ΜΑΧΑΙΡΙ Κ25, PHYTON, 32338

Designed in Spain
Tactical pocket knife
Brand: K25
Handle: Aluminium
Blade: Stainless steel (Camo coated)
Blade size: 13.5 cm
Blade thickness: 4 mm
Total size: 25 cm
Nylon sheath

Χαρακτηριστικά

The document was created on: Ιανουάριος 9, 2025